Getting started with the oximeter and AiLink

Read your measurements on the OLED screen, then connect AiLink if you want to view them on your phone.

Prepare the device

• Insert the two AAA batteries, following the polarity markings in the compartment.

• Read the precautions in your device’s manual before first use.

Take a spot measurement

• Insert a clean, dry finger. Turn on the oximeter and keep your hand still.

• Wait for the display to stabilise. Read the values directly on the screen; a phone is optional.

Identify the displayed values

• SpO_2: estimated oxygen saturation. PR: pulse rate. PI: perfusion index, when shown.

• This guide explains the controls and display, not the medical interpretation of a reading. The device is intended for spot measurements, not for recording an entire night.

AiLink · Guides

Connect AiLink

• Install “AiLink” from the App Store or Google Play on a compatible phone.

• Enable Bluetooth and grant the nearby-device permissions requested by the app.

• Turn on the oximeter near the phone. In AiLink, choose “Add Device”, then “OXIMETER” when it appears.

• Take a measurement and view the supported data in the app. Menus may vary by version.

If the connection fails

• Check that the oximeter is on, Bluetooth is enabled and AiLink has the required permissions. Restart the search in the app.

• The oximeter screen works without a phone. For support, provide your phone model, iOS or Android version and the step that fails. Hide personal readings in screenshots.

Keep your documents

• Keep the full manual and your unit’s packaging. Contact: contact@farmaqo.fr.

This quick-start guide supplements the model’s safety manual. Check the reference on your device and follow the instructions supplied with your unit.

FARMAQO practical guide · 29 September 2026 edition. Read alongside the manual supplied in the box. App screens may change.