Getting started with the thermometer and VitaConnect

Choose the appropriate mode, adjust the sound and view your measurements in VitaConnect.

Measure at the forehead

• Turn on the thermometer using the measurement button. Select adult or child forehead mode with the mode button.

• Aim at the centre of the forehead, between the eyebrows, from 0 to 3 cm away. Press the measurement button.

• Wait about 10 seconds between readings to avoid warming the sensor.

Switch to ear mode

• Remove the probe cap: ear mode is selected automatically. Position the probe without forcing it, following the manual, then press the measurement button.

• Follow the limits in the full FARMAQO manual: ear mode from one year of age; do not use on newborns.

Use object mode

• Select object mode, identified by the house icon. Aim at the surface from 0 to 3 cm away and press the measurement button.

• The illustrated guide shows a feeding bottle. This is a surface reading: never immerse the probe in bath water or another liquid. Do not use this mode for body temperature.

VitaConnect · Guides

Adjust sound and units

• With the thermometer on, touch the sound button to enable or mute the beep. The corresponding icon indicates silent mode.

• Hold the same button for 5 seconds to switch between °C and °F.

Connect VitaConnect

• Enable an Internet connection and Bluetooth on your phone. Install “VitaConnect” from the App Store or Google Play.

• Open the app, select sign-up and create an account with your email address and a password.

• Choose “Add device”, then the infrared thermometer. Turn on the thermometer and start the connection in the app. Wait for confirmation.

• Take a measurement using the device. The value then appears in the app. The supplementary leaflet also describes the history and profile under “Me”, then “Personal data”.

If you need help

• Check Bluetooth, app permissions and that the thermometer is on near the phone. Menus may change.

• Contact contact@farmaqo.fr with the thermometer reference, phone model and the step that fails. This guide supplements the manual’s precautions; it does not replace medical advice.

This quick-start guide supplements the model’s safety manual. Check the reference on your device and follow the instructions supplied with your unit.

FARMAQO practical guide · 29 September 2026 edition. Read alongside the manual supplied in the box. App screens may change.